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Common errors and troubleshooting immunoanalyzer (cobas)

Dr mehrdad vanaki

Consultant of QMS & QA in medical lab

Cobas e 411 immunoanalyzer ECL

cobas e 411 analyzer

Specifications

System	Fully automated, random access system for immunoassay analysis. It is available as both, a disk system and a rack system.
System components	Analytical module including Window XP embedded operated touch screen PC Sample handling module: rack or disk operated
Sample throughput	Up to 88 samples/hr (theoretical max)
Test throughput	Up to 88 tests/hr (theoretical max)
Number of channels	18 channels/reagent slots for up to 18 different assays
Programmable parameters	Max 60 assays definable via 2D-barcode (programming by loading)
Sample types	Serum, Plasma, Urine

Sample input/output	Load/unload capacity 30 samples (disk) 75 samples on 15 racks Rack RD standard 5 position rack Rack types Routine, STAT, Control, Calibrator STAT handling Any unoccupied position on the sample disk, dedicated STAT port on rack feeder
Sample container types	Primary tubes 5–10 ml; 16x100,16x75, 13x100, 13x75mm Sample cup 2,5 ml Micro cup not allowed Cups on tube Cup on top of a 16 x 75/100 mm tube
Sample volume	10 to 50 µl per test, depending on assay protocol
Minimum sample volume	Primary tubes : 600 µl (13mm tube), 1,000 µl (16mm tube) Sample cup: 200 µl (standard Cup on Tube) 150 µl with special settin
Sample data base	2.000 tests for routine, STAT and control results
Test methods	Pre-defined assay protocols (sandwich, competitive, titration)

ECL immunoassay (ECLIA) is based on the use of **ruthenium-tris(bypyridyl)** complex and **tripropylamine (TPA)** and the final chemiluminescent product is formed during detection step.


Figure B-2

Detection of a ruthenium-labeled immune complex


Figure B-3

The ECL reaction at the electrode surface

Three test principles are available on **COBAS e 411** analyzer:

a) Competitive principle:

- For extremely small analyte (e.g: FT3)

b) Sandwich principle:

- For larger analyte (e.g: TSH)

c) Bridging principle:

- For detecting antibodies in the sample



Figure C-1

ECL assay principles

Sandwich Principle (HBsAg)

FIRST REACTION

SECOND REACTION

LIGHT REACTION

Antigen

Ruthenium-labeled antibody

Biotinylated antibody

Streptavidin-coated microbead

TPA Tripropylamine

Pt. sample combine with reagent containing **biotinylated anti-HBsAg antibodies** and mixture of monoclonal/polyclonal anti-HBsAg labeled with ruthenium complex.

Streptavidin-coated microparticles added → the complex bound via interaction of biotin and streptavidin.

The rx mixture aspirated into measuring cell the microparticle magnetically bound on the electrode surface and unbound reagent washed away by ProCell


Application of voltage to the electrode stimulate the emission of light then the intensity of light emitted measured by photomultiplier. The light emitted **directly proportional** to the amount of HBsAg in the sample

Competitive principle (FT3)

FIRST REACTION

SECOND REACTION

LIGHT REACTION

TPA Tripropylamine


Sample and specific anti-T3 Ab labeled with ruthenium complex

Biotinylated-T3 and streptavidin coated magnetic microbead added → form antibody-hapten cplx → entire cplx bound to microbead thru rx of biotin and streptavidin

Rx mixture containing immune complex transported to measuring cell → those immune cplx are magnetically captured on electrode → unbound reagent + sample washed away by Procell

ECL reaction occurs when stimulated electrically, and the light produced is **indirectly proportional** with [] of Ag in pt. sample.

Bridging Principle (IgG/IgM)


Similar with sandwich principle except the assay designed to detect antibodies, not antigens. This accomplished by including biotinylated and ruthenium-labeled antigens in the reagent for which Ab has affinity.

Example

1 step test (sandwich assay)

Principle: one-step sandwich assay

Incubation (9 minutes)

10 μ L of sample, a biotinylated monoclonal hCT-specific antibody and a monoclonal hCT-specific antibody labeled with ruthenium complex are incubated. A sandwich complex is formed with hCT carrying a biotinylated and a ruthenylated antibody.

Incubation (9 minutes)

After addition of streptavidin-coated microparticles, the complex becomes bound to the solid phase via interaction of biotin and streptavidin.

Measurement

The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.

cancer.

Test principle: one-step sandwich assay

Biotinylated monoclonal antibody against human HE4

Ruthenylated monoclonal antibody against human HE4

Streptavidin microparticle

Elecsys technology

ECL (ElectroChemiLuminescence) is Roche's technology for immunoassay detection. Based on this technology and combined with well-designed, specific and sensitive immunoassays, Elecsys delivers reliable results. The development of ECL immunoassays is based on the use of a ruthenium-complex and tripropylamine (TPA). The chemiluminescence reaction for the detection of the reaction complex is initiated by applying a voltage to the sample solution resulting in a precisely controlled reaction. ECL technology can accommodate many immunoassay principles while providing superior performance.

Reagents

One grip handling of reagents
3 ways to unmatched performance

As the heart of the cobas e 411 analyzer, ECL technology delivers the unique combination of nine minute STAT applications, high analytical sensitivity and wide measuring ranges for unparalleled performance in

2- High analytical sensitivity

enables the use of innovative tests and requires less sample volume

examples:

- Troponin T 4th gen. with $<10\%$ CV at 95th population percentile
- hCG+, and TSH requires 10 and 50 μL sample volume, respectively

1-Wide measuring ranges reduces the amount of tests rerun due to out of range results

examples:

- Estradiol range of 5 – 4,300 pg/mL
- Progesterone range of 0 – 60 ng/mL
- total-PSA range of 0 – 100 ng/mL

s

Features


ECL technology

Virtually eliminates repeats and reruns :

The measuring cell uses an amplified signal to detect ultra-low concentrations of analyte. By precisely controlling electrically initiated reactions, ECL technology virtually eliminates unnecessary repeats and reruns, providing accurate results the first time

"Electro" refers to electrical stimulation.
+
"Chem" indicates a chemical reaction.
+
"Luminescence" means "produces light."
=
Electrochemiluminescence (ECL)

Reaction Phase – Measuring Cell

SYSTEM REAGENT

Procell

- Conditioning of the Electrode.
- Transport of the assay reaction mixture
- Washing of the streptavidin coated microbeads

Cleancell

- Cleaning of the tubing system and of the measuring cell after every measurement.

Conditioning of the Electrode.

Sysclean

Sodium Hypochlorite Solution used for the cleaning of the measuring cell
(every two weeks)

Syswash

Increasing rinsing efficiency between pipetting steps

Preventing Bacterial Growth

Syswash added to the Distilled water container, with a dilution of 1+100.

System Reagent

Procell :

- 1- conditioning of the electrode
- 2-Transport of the assay reaction mixture
- 3- Washing of the streptavidin microbeads

Cleancell

Cleaning of the tubing system and the measuring cell after every measurement

Table of contents

Revision history	2
Contact addresses	3
Table of contents	5
Preface	7
Where to find information	7
Online Help system	8
Customer information	9
Conventions used in this manual	9
System description	Part A
1 Safety information	
Safety classifications	A-5
Safety precautions	A-5
System safety labels	A-18
2 System overview	
Overview of the cobas e 411 analyzer	A-27
System components	A-30
Analyzer unit	A-30
Control unit	A-34
Power components	A-34
System specifications	A-37
3 Control unit	
Components	A-55
Software basics	A-60
4 Analyzer components	
Overview of the analyzer components	A-71
Sample/reagent area components	A-72
Consumables area components	A-89
Measuring area components	A-97
Operation	Part B
5 Overview	
Operation Guide	B-5
Disk and rack systems	B-5
Software keyboard	B-5
6 Daily operation	
Overview	B-11
Pre-start inspection	B-12
Pre-routine operation	B-18
Routine operation	B-31
Results	B-81
Post-operation data management	B-87
Daily maintenance	B-87
Switch off the analyzer	B-88
7 Special operation: How to...	
Overview	B-93
How to apply barcode labels to tubes and vials	B-94
How to upload results manually	B-96
How to select two or more samples in Data Review	B-98
How to archive sample data	B-99
How to rerun a sample	B-100
How to reset the system	B-101
How to delete a single open request	B-103
How to assign calibrator and control positions	B-104
How to install and define Roche controls	B-106
How to install and define non-Roche controls	B-107
How to edit control values	B-109
How to request standby bottle QC	B-110
How to cancel a print job	B-111
How to print a daily Alarm Trace report	B-112
How to print a System Communication Trace report	B-113
How to change documentation settings	B-114
How to program calculated tests	B-115
How to change test settings	B-116
How to change calibration settings	B-118
How to change the sample disk mode	B-119
Maintenance	Part C
8 Maintenance	
Overview	C-5
Maintenance schedule	C-8
Maintenance log	C-9
Daily maintenance	C-12
Weekly maintenance	C-15
Every two weeks maintenance	C-19
Every two months maintenance	C-24
As needed maintenance	C-30
Troubleshooting	Part D
9 Troubleshooting	
Introduction	D-5
Alarms	D-5
Problem categories	D-5
Contacting Technical Support	D-8
Basic troubleshooting flow chart	D-10
Immunoassay troubleshooting	D-11
Instrument troubleshooting	D-14
Chemistry troubleshooting	D-24

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Sample/reagent problems

- Insufficient sample volumes of:
 - Sample cups
 - Sample tubes
-  See:
 - Residual volume (disk system)* on page A-42
 - Residual volume (rack system)* on page A-42
- Reagent packs not loaded correctly with the white bottle furthest from the center of the reagent rotor
- Diluent bottles not loaded with the correct barcode label facing outwards
- Incorrect sample preparation
 - Proper reconstitution
 - Aliquoting/labeling
- Incorrect sample storage
 - Refrigerating
 - Freezing
 - Thawing
 - Labeling after use
- Manually program calibrators
- Operator error

Instrument/hardware problems

- Electrical/electronic
- Mechanical:
 - Pinch valve tubing leaking
 - AssayCup/AssayTip jammed
 - Problem with solid waste tray removal
- Operator error:
 - Incorrect size of sample tubes or cups used
 - Incorrect replacement of AssayCup/AssayTip
 - Incorrect or insufficient maintenance
 - Caps not removed from liquid waste containers
 - Caps not opened on ProCell/CleanCell bottles

Computer/software problems

- Incorrect system parameters
- Faulty barcode card or software settings
- Disk system:
 - Single disk mode
 - Multiple disk mode
 - Starting from Standby compared with starting from S. Stop

Facility problems

- Temperature
- Humidity
- Power supply
- Water supply
- Drain

👁 See *System specifications* on page A-37.

Operator's primary responsibility

Your primary responsibility lies in the following areas:

- Reagent, calibrator, and control preparation and storage
- Sample and sample container preparation
- Computer parameters and general computer input/output operations
- Basic operator technique, including instrument and computer operation
- Basic component replacement
- Maintenance

For essential safety information, see Electrical safety precautions on page A-7.

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