# VIDAS®

- **VIDAS®** is a reliable and easy-to-use automated benchtop immunoanalyzer. Based on the Enzyme Linked Fluorescent Assay (**ELFA**) technology, it provides high quality on-demand test results.

- **Principle of a VIDAS® Assay**

The Solid Phase Receptacle (SPR®) serves as the solid phase as well as the pipetting device for the assay. Reagents for the assay are ready-to-use and predispensed in the sealed reagent strips. All of the assay steps are performed automatically by the instrument. The reaction medium is cycled in and out of the SPR® several times.

Unbound components are eliminated during the washing steps. During the final detection step,

the substrate (4-Methyl-umbelliferyl phosphate) is cycled in and out of the SPR®. The conjugate enzyme catalyzes the hydrolysis of this substrate into a fluorescent product (4-Methyl-umbelliferone), the fluorescence of which is measured at 450 nm.

The intensity of the fluorescence depends on the concentration of alkaline phosphatase present on the SPR® that transforms the substrate. At the end of the assay, results are automatically calculated by the instrument. For some tests, two detection steps are performed successively.

For antigen detection, the SPR® is generally coated on the interior with capture antibody or sometimes with a derivative of the analyte. For antibody detection, the SPR® is coated with a capture antigen or antibody directed to the antigen. Depending on the test, the conjugate can be a derivative of the analyte or an antibody labeled with alkaline phosphatase. For more details, refer to the assay package inserts

# VIDAS®

# 3

## System Description and Basic Operations

### *mini VIDAS®* analyzer Description

Figure 3-1: *mini VIDAS®* analyzer

- 1 — Operating status light
- 2 — SPR® block
- 3 — Strips section
- 4 — Strip preparation tray
- 5 — Keypad and screen
- 6 — Thermal printer
- 7 — External bar code reader

**Note:** Previous versions of the *mini VIDAS®* analyzer may not have a strip preparation tray.

# VIDAS® *Because precision helps*

## *VIDAS® instruments, a Customized solution for your daily workload*

VIDAS® instruments are bench-top automated immunoassay systems based on ELFA\* technology. All VIDAS® instruments are operational 24/7 and use the same reagent references.

They allow you to process several different parameters simultaneously, in both single random access sample and batch test mode, for all types of analyses. Their unique design guarantees that there is no inter-reagent or inter-sample contamination.

**miniVIDAS®**  
12 tests simultaneously  
Integrated software and printer  
Limited space required due to its compact size  
Load and Go

**VIDAS®**  
30 tests simultaneously  
Ergonomical and intuitive software

*Robust instrument  
No daily maintenance  
Limited training time*

\* Enzyme-linked fluorescence assay.

## EASY-TO-USE, OPERATIONAL 24/7, JUST LOAD AND GO

Patient ID acquisition & sample distribution

Place SPR® and STRIP into the machine


Press Start key

# 1- Measuring Principle

The reaction occurs within the interior of the SPR whereby antibodies and conjugate form a sandwich.

(4-MUP) is cycled into SPR and conjugate enzyme catalyses the hydrolysis of the substrate into 4-Methyl-umbelliferone which is measured at 450nm.

The figure below represents a typical label:


Figure 3-7: Reagent strip label

- 1 — Alphanumeric bar code.  
The first two characters of the alphanumeric code indicate the unique assay code (UAC). If the system is unable to read this code, the alphanumeric form of the bar code is used to enter the bar code.
- 2 — Assay code  
The assay code is based on three or four characters which represent the name of the assay to be performed.  
Examples:  
RBG = RUB IgG (Immunoassay)  
153 = CA 15-3  
All references to an assay in the software and on printed reports use the assay code.
- 3 — Lot number  
The lot number is used by the **mini VIDAS®** analyzer to code the expiration date of the reagent lot.
- 4 — Expiration date  
This corresponds to the expiration date of the reagent lot.  
When the bar code reader reads this date, the **mini VIDAS®** analyzer compares it to the date stored in its memory to ensure that the lot has not passed its expiration date.
- 5 — Sequence number  
Within a single batch, each reagent tray carries a different number corresponding to the order of manufacture.
- 6 — Bar code  
Read by the internal bar code reader.
- 7 — ID  
Optional: A blank space is provided to write the sample ID (if required).

The SPR® is a plastic (polystyrene) device capable of capturing soluble proteins, viruses and bacteria. It is sealed with a color-coded, barcoded dot perforated in the center. Each SPR® is disposable.

It eliminates cross-contamination between reagent and instrument and reduces maintenance to a minimum since there is no tubing, syringe or sampling needle involved.

Each SPR® for a specific assay is identified by a color-coded dot on which a bar code and a letter and/or digit code representing the assay name are printed.

**VIDAS®** clinical assays

**VIDAS®** industry assays

Figure 3-8: Solid phase receptacle (SPR®)

The SPR® is the solid phase base for the immunological reaction. Its interior walls are coated with antibody or antigen that captures a target analyte.

The target analyte from the sample binds to the SPR®'s interior coating (antibody, antigen, etc.). It is then bound by an enzyme-conjugated antibody or antigen, forming a «sandwich». The immobilized enzyme catalyzes the hydrolysis of the substrate into a fluorescent end product.

- 
- The SPR<sup>®</sup> is used to pipette samples and reagents and perform the following operations:
    - sampling
    - incubation
    - mixing
    - washing

\* SPR : Solid Phase Receptade

SPR\*

EMISSION

SUBSTRATE

WASHING BUFFER

RECEPTION


DILUENT

WASHING BUFFER

DILUENT

SAMPLE


REAGENT STRIP


# VIDAS Principle

## CAPTURE OF ANTIGENS


The antibody captures the target pathogens

## SANDWICH TEST



A second antibody conjugated with an enzyme binds to specific antigens

## DETECTION

The intensity of the reaction is interpreted by the system

# VIDAS LMX PRINCIPLE

Specificity

+

Sensitivity

=

Next Day Results

## Analysis Methods

The *mini VIDAS*<sup>®</sup> analyzer uses several different methods to calculate results. Three basic categories of analysis methods are:

Calculation method	Description
Test sample to Standard	A comparison of the Relative Fluorescence Value (RFV) of the test sample to that of a standard. This method is used in Single Reagent Strip qualitative and most semi-quantitative assays.
Test sample to Reference	A comparison of the RFV of the test sample to that of a reference. This method is used in Dual Reagent Strip qualitative assays.
Curve Fitting Equations	The RFV of a test sample is mathematically placed on a Calibration Curve. This method is used for all quantitative and some semi-quantitative assays.

## ***Qualitative Assays (Single Reagent Strips)***

Two analysis methods are used for single reagent strip qualitative assays: the P/S method and the P-S method.

In these methods, the “P” stands for the RFV of the test sample and the “S” stands for the RFV of a standard.

### **Assay Calibration**

Single reagent strip assays are calibrated using one or two standards/calibrators supplied with the assay kit. Calibration using the standard(s)/calibrator(s) provided in the kit should be run the first time the assay kit lot is used, after the master lot data have been entered.

If replicate standards are run, their values are averaged. The calibration obtained can be used for a programmed period of days.

After that, the software automatically expires the calibration, requiring that you run another one from the same lot.

## P/S Method

For assays that use the P/S method, the test value is calculated from the ratio of the sample RFV to that of the standard RFV.

- RFV (tested sample) = 2158
- RFV (standard) = 2177
- Test value (VT) =  $2158/2177 = 0.99$

## P-S Method

For assays that use the P-S method, the test value is calculated from the difference between the sample RFV and the standard RFV.

- RFV (tested sample) = 1774
- RFV (standard) = 1689
- Test value (VT) =  $1774 - 1689 = 85$

## Assay Result

In both the P/S and P-S methods, the result of the assay is determined by comparing the test value to a set of preprogrammed thresholds.

Depending on the assay and the test value result compared to the thresholds, the result can be either positive, equivocal or negative.

For every assay type, there is a high and low threshold value. The assay result is interpreted from the test value as follows:

If the test value is:	the result is:
$\geq$ high threshold	Positive
$<$ high threshold and $\geq$ low threshold	Equivocal
$<$ low threshold	Negative

**Note:** For some qualitative assays (e.g. TOXO Competition, anti-HBc total), the result of the test value will be negative if the test value is to the high threshold.

# Quantitative Assays

**Quantitative VIDAS® assays use a calibration curve to determine analyte concentrations in test samples. Note:**

**Although standards/calibrators can be run after analyses have completed, it is recommended that they be run either before or with analyses from the same assay kit lot. If an assay has two standards/calibrators, both must be run at the same time**

## Principles Of Calibration / Recalibration

The reagents are calibrated according to two different modes, depending on whether the assays are quantitative or qualitative.

Assay	Description
Quantitative	Calibration is performed at the factory for each new reagent lot and is entered into the instrument using a bar code provided with each reagent kit. This creates a «Master Curve» that is stored in memory. It can be readjusted by running a standard/calibrator. Readjustment of the curve must be validated by testing the control(s) in the kit. Recalibration should be performed every 14 days or 28 days depending on reagents.
Qualitative	The result is interpreted by comparing the «Test Value» to one or two threshold values entered into the <i>mini VIDAS</i> ® software.

### Calibration

#### Theoretical Principle for Quantitative Assays

## Calibration

### Theoretical Principle for Quantitative Assays

The principle consists of determining the mathematical equation which represents the calibration curve, i.e. the relationship between the RFV and the concentration of the standards/calibrators.

Well-defined reference solutions (standards) are used to make this determination. This calibration curve is established with at least 5 standards/calibrators. It is valid between zero and the standard/calibrator with the highest titer.

3 mathematical models can be used to establish the master lot curve:

1. 4 logistic parameters or Rodbard model

$$\text{RFV} = \frac{a_1 - a_4}{1 + (c / a_3)^{a_2}} + a_4$$

2. Polynomial

$$\ln(c) = a_1 + a_2 \text{ RFV} + a_3 \text{ RFV}^2 + a_4 \text{ RFV}^3$$

3. Semi-log

$$\text{RFV} = a_1 + a_2 \ln(c)$$

c represents the concentration.

$a_1$ ,  $a_2$ ,  $a_3$ ,  $a_4$  are the mathematical parameters of the model. They are recalculated for each calibration.

## Principle Applied to the *mini VIDAS*® analyzer

For the *mini VIDAS*® analyzer, calibration is performed during the production of each new lot of SPR®s and reagents. The number of standards/calibrators varies between 5 and 11 depending on the biological assay. Each lot is associated with a particular mathematical model.

To determine the master curve, the standards/calibrators are tested in seven different runs on the same *mini VIDAS*® analyzer. The mean curve of these seven runs becomes the master curve.


Figure 3-10: Master curve (example)

The master curve is specific for a certain manufacturing lot and biological parameter.

## Principle of Use

Upon receipt of a new lot of reagents, the user must enter the master curve (master lot data) using the Master Lot Data (bar code) provided with the kit. The standard/calibrator should then be run to readjust the curve. The controls are tested in the same run to check the recalibration.

Recalibration is valid for 14 days or 28 days depending on reagents, after which, the software requests another calibration.

To avoid any deviation, the master lot curve is the one set in the factory and not the last curve readjusted. The software can perfectly manage simultaneous use of several master curves corresponding to different lots for a given assay.


Figure 3-12: Master curves

## Establishing the Master Curve

Each kit includes the master curve in the form of:

- an MLE bar code printed on the reagent kit label,
- an MLE card provided in the reagent kit (for certain assays).

The bar code contains the following information:

- 4 RFV levels, corresponding to 4 standards defined in the software. The master curve is restored using 4 points. The levels of concentration of these 4 standards are known to the **mini VIDAS®** software and are specific for each assay.

The **mini VIDAS®** software replots the total curve using the 4 pairs of concentrations and RFVs (master lot data).

The remaining information includes:

- the assay code,
- the kit lot number (SPR®s and reagents). It is used to link the master curve data to the strip lot number,
- mathematical model used to establish the master curve (7: Rodbard, 8: Polynomial, 9: Semilog),
- the dose value for the standard (recalibrator) which is essential for calculating the correction factor. This value may vary slightly from one lot to another,
- range values for the controls in the kit,
- range values for the standard's RFV,
- maximum variation coefficient of the standard's RFV duplicates or triplicates.

All these data are specific to an assay and a given lot.


The standard/calibrator must be identified by S1 or S2 (e.g. HBET) or S1 and S2 (e.g. TSH3).

Calibration should be performed each time a new lot is opened (after master lot data have been entered), and every 14 or 28 days, depending on the reagent.


## Initial biological investigation of thyroid disorders

### Patient not receiving treatment


## Biological monitoring\*

### Patients treated with thyroxin as a substitute


\* A considerable time period (at least 2 to 3 weeks) should separate the

### Patients treated with synthetic anti-thyroid drugs



biological follow-up from the initiation or modification of treatment.


## Two-step indirect sandwich assay format

- Monoclonal antibodies anti-p24
- Synthesized peptides to capture HIV-1, 2+0 antibodies
- Biotinylated p24 polyclonal antibodies
- Biotinylated antigens
- p24 Ag in the sample
- Patient antibodies in the sample

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---

Maintenance of the **mini VIDAS<sup>®</sup>** analyzer comprises:

- the preventive maintenance operations performed by bioMérieux technicians or a qualified person trained by bioMérieux,
- the routine maintenance operations which should be performed by the user.

### Preventive Maintenance

---

***IMPORTANT: Optimum performance of the **mini VIDAS<sup>®</sup>** analyzer (within the limits of bioMérieux specifications) depends on the preventive maintenance being performed by bioMérieux technicians or a qualified person trained by bioMérieux only, as part of a maintenance contract.***

# A

## Appendix - Maintenance Records

---

### Maintenance Schedule

---

Part	Frequency	Procedure
SPR® block	Monthly	Cleaning the SPR® block (Monthly)
Reagent strip tray	Every 6 months or if needed	Cleaning the Strip Preparation Trays (Every 6 Months or if Needed)
Plastic tray (underneath the reagent strip tray)	Every 6 months or if needed	Cleaning the Plastic Trays (Every 6 Months or if Needed)
Strip preparation tray	Every 6 months or if needed	Cleaning the Strip Preparation Trays (Every 6 Months or if Needed)
Outer covers	Every 6 months or if needed	Cleaning the Outer Covers of the Instrument (Every 6 Months or if Needed)
Screen and keypad	Every 6 months or if needed	Cleaning the Screen and Keypad (Every 6 Months or if Needed)
Bar code reader	Every 6 months or if needed	Cleaning the Bar Code Reader (Every 6 Months or if Needed)
Optical lenses	Monthly	Cleaning Optical Lenses (Monthly)
Pipetting device	Monthly	Checking the Pipetting System (Monthly)

## WEAR

[illegible]


## Cleaning the SPR<sup>®</sup> block (Monthly)

1. Shut down the *mini VIDAS*<sup>®</sup> analyzer (see [Shutting Down the System on page 4-4](#)).

**Note:** Previous versions of the **mini VIDAS®** analyzer have no safety latch. For these versions, you can directly open the SPR® block door

4. Using a Dacron<sup>®</sup> swab moistened with a laboratory cleaning and disinfection solution, carefully clean the interior of each SPR<sup>®</sup> liner of the first SPR<sup>®</sup> block.
5. Continue the operation by using a new Dacron<sup>®</sup> swab moistened with the cleaning and disinfection solution for each SPR<sup>®</sup> block.

**Note:** *Laboratory cleaning and disinfection solution can be identical to the one used in the wipes (same active principle).*

**IMPORTANT:** *The presence of fibers from the Dacron<sup>®</sup> swab inside the SPR<sup>®</sup> liners may interfere with the performance of the mini VIDAS<sup>®</sup> analyzer.*

**When cleaning the SPR<sup>®</sup> liners, replace any swab if it becomes fluffy.**

6. Clean all surfaces of the SPR<sup>®</sup> block by using the cleaning and disinfection wipes.

**Note:** Wipe all surfaces thoroughly with the wipe, ensure complete wetting and allow to take effect.

**CAUTION:** Do not clean the seals situated above the SPR<sup>®</sup> block as they only require cleaning when preventive maintenance is performed by the bioMérieux Field Service Engineer.


In the event of accidental spills or contamination, it is imperative that you contact bioMérieux or your bioMérieux representative.

Handling of the seals may affect the performance of your instrument.

---

7. Press each SPR<sup>®</sup> liner until you are able to clean its rear end.
8. Repeat this procedure for each SPR<sup>®</sup> liner of each block.
9. If no other maintenance tasks have to be performed, start the **mini VIDAS<sup>®</sup>** analyzer (see [Starting the System on page 4-3](#)).

## Cleaning Optical Lenses (Monthly)

It is recommended to perform the following procedure to ensure the cleanliness of the scanner and prevent optical self-test errors and self-calibration drifts.

**IMPORTANT:** *The VIDAS<sup>®</sup> Lens Cleaner is intended for cleaning the optical lenses of the mini VIDAS<sup>®</sup> analyzers.*

- *Turn the power switch to OFF before cleaning the optical lenses.*
- *Before each use: Take the VIDAS<sup>®</sup> Lens Cleaner out of its package and squeeze the rubber ball hard two or three times in order to remove any moisture inside.*
- *During cleaning: Avoid touching the lenses with the extremity of the tool as it could damage them.*
- *After each use: Put the VIDAS<sup>®</sup> Lens Cleaner back into its package and store it in a dry place.*
- *In case of accidental spills on the VIDAS<sup>®</sup> Lens Cleaner, clean the tool with a dry clean cloth.*
- *If any substance is accidentally introduced inside the VIDAS<sup>®</sup> Lens Cleaner, replace the tool.*

1. Open the SPR<sup>®</sup> block doors of sections A and B.

The Scanner head should be visible (see figures below) either in section A or section B.

There are two scanner head models (see [Figure 8-1](#)).

- ***If any substance is accidentally introduced inside the VIDAS® Lens Cleaner, replace the tool.***


1. Open the SPR® block doors of sections A and B.

The Scanner head should be visible (see figures below) either in section A or section B.

There are two scanner head models (see [Figure 8-1](#)).

**Model A**

**Model B**

Figure 8-1: Scanner models

2. Shut down the ***mini VIDAS®*** analyzer (see [Shutting Down the System on page 4-4.](#)).

3. Position the **VIDAS®** Lens Cleaner above the lens.
4. Clean the upper lens (*tilted surface*): squeeze the rubber ball hard ten times.

**Model A**

**Model B**

Figure 8-2: Cleaning the upper lens

5. Clean the lower lens (*horizontal surface*) by squeezing the rubber ball ten times.

**Model A**


**Model B**

Figure 8-3: Cleaning the lower lens.

### ***Cleaning the Reagent Strip Trays (Every 6 Months or if Needed)***


1. Shut down the **mini VIDAS®** analyzer (see [Shutting Down the System on page 4-4.](#)).
2. Open the section tray dust cover (Plexiglas).

### ***Cleaning the Plastic Trays (Every 6 Months or if Needed)***

1. Shut down the **mini VIDAS®** analyzer (see [Shutting Down the System on page 4-4.](#)).
2. Open the section tray dust cover (Plexiglas).
3. Move the strip tray in order to have access to the plastic tray.

### ***Cleaning the Outer Covers of the Instrument (Every 6 Months or if Needed)***

1. Shut down the **mini VIDAS®** analyzer (see [Shutting Down the System on page 4-4.](#)).
2. Wipe all surfaces with cleaning and disinfection wipes.

**Note:** *Wipe all surfaces thoroughly with the wipe, ensure complete wetting and allow to take effect.*

3. If no other maintenance tasks have to be performed, start the **mini VIDAS®** analyzer (see [Starting the System on page 4-3.](#)).

### ***Cleaning the Screen and Keypad (Every 6 Months or if Needed)***

1. Shut down the **mini VIDAS®** analyzer (see [Shutting Down the System on page 4-4.](#)).
2. Wipe the screen and the keypad with cleaning and disinfection wipes.

**Note:** *Wipe all surfaces thoroughly with the wipe, ensure complete wetting and allow to take effect.*

3. If no other maintenance tasks have to be performed, start the **mini VIDAS®** analyzer (see [Starting the System on page 4-3.](#)).

### ***Cleaning the Bar Code Reader (Every 6 Months or if Needed)***

1. Wipe the bar code reader (except the window) with cleaning and disinfection wipes.

# ADVANCES IN IMMUNOASSAY TECHNOLOGY

Copyright © 2012 InTech

**INTeCH**

open science | open minds

free online editions of InTech  
Books and Journals can be found at  
[www.intechopen.com](http://www.intechopen.com)

## Part 1 New Materials

Chapter 1 Recombinant Antibodies  
and Non-Antibody Scaffolds for Immunoassays 3  
Bhupal Ban and Diane A. Blake

Chapter 2 Polyacrylonitrile Fiber as Matrix for Immunodiagnosics 23  
Swati Jain, Sruti Chattopadhyay,

## Part 2 Label-Free Technologies

Chapter 4 Fundamentals and Applications of Immunosensors 65  
Carlos Moina and Gabriel Ybarra




Fig. 2. Scheme of the basic integrated units that conform a biosensor.

## Part 3 Multiplexing Technologies

# Luminex:

## xMAP (multiplex analyte profiling) technology

# Luminex:

## xMAP (multiplex analyte profiling) technology