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standard process of abnormal results in medical labs

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suggested comment: **rechecked with records of patient /rechecked**

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suggested comment: **rechecked with same sample & new method (same method) /
please repeat with new sample**

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suggested comment: **rechecked & confirmed with same sample & new sample & new
method (same method) /**

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example:HIV Ab is suspicious with screening test with two samle & two method
(ELISA method & CLIA) & confirmed result by western blot method (as a
confirmatory test)

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Abnormal results / infectious panel – cardiac & tumor marker pannel

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- This sample aged(not fresh) sample / please repeat test with new fresh sample
- Collecting water mouth instead of a sputum sample/ final results of sputum smear & culture is not reliable & Interpretab
- Sample was milky / turbid / semiturbid /hemolysis (slightly – moderate- severe) / results is not reliable & Interpretable/ specialy for this test lactate/ ammonia / K / LFT & cardiac panel
- Sample was milky / turbid / semiturbid / results of lipid profile is not reliable & Interpretable
- Failure to the cold chain for carrying samples so final results does not have interpretive value

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suggestion of confirmatory tests for screening test

- this method is screening method with high sensitivity & low specificity with high false positive results
- RF latex method (serologic & screening method with 40% false positive rate) is not specific tests for RA disease & is often seen in cases of chronic infections & other systemic inflammatory conditions so Anti MCV & Anti-CCP & RF(IgG & IgM) with high specificity as confirmatory tests is highly recommended
- The positive titre of wright or widal Test is not a sign of acute illness so repeat of titer assay after four week is highly recommended / for more & fast evaluation of infection brucella (IgG & IgM) as confirmatory test is recommended.
- False positive results of PSA :
- Manipulation & massage prostate before examination
- Severe inflammatory disease & viral disease(UTI / prostatitis/ influenza ...)
- Intercourse with partner 24 before examination

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suggestion of confirmatory tests for screening test

- suspicious results of HIV Ab(with two new sample & new method) by elisa screening method is not reliable & interpretable so
- HIV PCR & HIV western blot as confirmatory tests is highly recommended
- suspicious (reactive) results of HCV Ab(with two new sample & new method) by elisa screening method is not reliable & interpretable so
- HCV PCR & HCV RIBA as confirmatory tests is highly recommended

شرح	کد متن آماده
Polynucleosis (Lekocytosis) observed.	۱
Anti Nuclear Ab is Positive , titer is : 1/20 (Normal Range : Up to 1/20) ANA Pattern : Homogenous Pattern .	۲
Urea Breath Test (UBT) for monitoring of Antibiothraphy is recommended .	۳
(Ferritin) With Same Sample & Chemiluminescence Method Rechecked & Confirmed .	۴
Two Time Whith Two Sample Rechecked.	۵
30 - 40 % (False Negative & Positive for VDRL & RPR) in Collagen Vascular Disease & chronic inflamation observed. Please Repeat VDRL after four week.	۶
Confirmatory test , such as , PCR . RIBA ect . Is highly Recommended .	۷
Repeated sample of 14 days Interval is recommended .	۸
Using 3rd generation Elisa kit , the sample has a Positive reaction , for more evalution PCR is recommended .	۹
Pancytosis (Leukopenia & Anemia & Thrombocytopenia) Bipicture blood morphology (Macrocyte & Normochrom , Microcyte , Hypochrom) Lymphocytosis with Smudge Cell (Moderate) observed.	۱۰
Severe Lymphocytosis (Mature Lymph)	۱۱
Cytochemical Stain (MPO & SSB) Other confiromatory test (Bone Marrow ASP & Flowcytometry) is Highly recommended .	۱۲
After 50 gram oral glucose	۱۳
After 50 gram oral glucose . Abnormal GCT is above 140 mg/dl and GTT is required	۱۴
After 75 gram oral glucose . The higher the level of two or more than two measurments , favor the diagnosis of gestational diabetes.	۱۵
After 100 gram oral glucose . The higher the level of two or more than two measurments , favor the diagnosis of gestational diabetes.	۱۶
Serum Iron & Ferritin & Hb Elec.. is Highly recommended .	۱۷
Severe Oligospermia observed.	۲۶
This SCREENING test shown to have False Positive (7% to 10%) and Negative(2%) result in local and international studies . DIAGNOSTIC procedure to obtain conclusive result is Highly recommended. Risk factor in another calculation software (Prisca) shows Trisomy 21 risk for this	۲۹
LDL has done with Quantitative chemistry method .	۳۰
Azo Spermia (Dont See any Sperm after centrifugation) .	۳۲
(T4) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۳
(T3) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۴
(TSH) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۵
(T4 & TSH) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۶
(PSA) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۷
(Prolactin) With Same Sample & Electrochemiluminescence Method Rechecked & Confirmed .	۳۸
A Fresh Sample For More Evaluation is Recommended .	۳۹
Spermogram Two Time With New Sample Reched Oligospermia observed.	۴۰
Variant Lymphocyte (Atypical Lymph) : (<6% Normal) Variant Lymphocyte form is a new terminology for atypical lymphocyte . expected level for normal individual is <6%	۴۱
Hb Electerophoresis is Highly Recommended .	۴۲
	۴۳

Severe Shift to the Left (myeloid series) have seen (Bone marrow Asp)

Platelet Aggregation & Giant Platelet observed.

Hypersegmented Neutrophil observed.

Urine Culture is Highly Recommended .

Anti sperm Ab is highly Recommended .

Serum Was Milky .

Serum Was Semi Turbid .

Serum Was Slight Turbid .

Serum Was Turbid .

Lymphocyte Variant From (LVF) is a New Terminology For Atypical Lymphocyte .

Probably Contamination Repeat of Urine Culture is recommended .

Uncomputable for TG over >400 mg/ml

Scotch Tape Test Negative .

Rouleaux Formation observed.

Differentiative index of minor thalassemia & Iron Deficiency

MCV/RBC (Mentzer Formula) :

a.Iron Deficiency Anemia ≥ 13.0

b.Beta-Thalassemia Minor < 13.0

MCH/RBC (Srivastava & Benington) :

a.Iron Deficiency Anemia ≥ 4.4

b.Beta-Thalassemia Minor < 4.4

MCH/RBC (Srivastava & Benington) :

a.Iron Deficiency Anemia ≥ 4.4

b.Beta-Thalassemia Minor < 4.4

B12 & Folate Serum for Screening of Megaloblastic Anemia is Suggested .

The Software estimates prenatal screening risk , between the gestational age of risks for this patient is not possible .

OA (Ab) Probably is cross reactive Antibody.

Please repeat Widal test with new sample after (4) four week .

Probably past vaccination or postinfection .

The result must be confirmed by PCR or Western Blot Analysis if Clinically indicated .

The Test should be rechecked with afresh sample .

(T3 & T4 & TSH) With Same Sample & Chemiluminescence Method Rechecked & Confirmed .

Our lab are not responsible for patient identification .

Eosinophilia observed.

Serum IgE is Recommended .

(HBS Ag)With Same Sample & Chemiluminescence Method Rechecked & Confirmed

Toxic Granulation observed.

For recent infection , HAV IgM is recommended.

Dipicture RBC Morphology observed.

Repeated sample of 14 days interval is recommended .

ANA Ab were seen , titer is : 1/80

(Normal range : Up to 1/20)

Homogenous & Nucleolar pattern .

Kidney microsomal Ab were seen , titer is : 1/40

Anti Smooth muscle Ab is positive , titer is : 1/20

(Normal range : Up to 1/10)

Urine Culture before antibiotherapy is recommended .

GTT for confirmation of gestational diabet (GD) is recommended .

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Summary

- While many steps are involved in method validation, not all have to be performed by the laboratory
- Lab must verify accuracy, precision, AMR, CRR, and reference intervals before methods are adopted (and must be signed by lab director)


Fig. 2. Validation requirements according ISO 15189:2007.

- Non standard method
- Laboratory designed by developed method
- Standard method used outside their intended scope
- Modified validated method

VALIDATION

Before use as diagnostic test/method

DEFINE PERFORMANCE CHARACTERISTICS,
As extensive as is necessary to confirm,
Through the provision of objective evidence,
That the specific requirements for the intended use have been fulfilled

DOCUMENT and record the results obtained and the
procedure used for the validation

CONTINUOUS VALIDATION

- Existing METHOD with defined performance
- Existing METHOD used after repair

VERIFICATION

Before use as diagnostic test/method

COMPARE PERFORMANCE CHARACTERISTICS, with
specifications.
Be aware of and respect specified requirements

DOCUMENT and record the results obtained and the
procedure used for the validation

CONTINUOUS VALIDATION

VALIDATION VS VERIFICATION

OF A METHOD IN A CLINICAL/MEDICAL LABORATORY

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Verify Reportable Range

- Reportable range includes:
 - Analytical measurement range (**AMR**) – range of values that instrument can report directly (less accurately called **linearity**)
 - Clinically reportable range (**CRR**) – range of values that can be reported with dilution or concentration of samples
- Retain records while method in use, at least 2 years after discontinued

Verify Reportable Range

- **AMR** must be verified before method introduced, and at least every 6 months (and after recalibration or major maintenance) while in use
- **CRR** is a clinical decision by lab director, and does not require experiments or re-validation; however, dilution or concentration protocols must be specified in methods

CAP (two concepts) CLIA 493.2

Analytical measurement range: range of analyte values that a method can directly measure on the specimen without any dilution, or other pretreatment not part of the usual assay process.

Reportable range: span of test result values over which the lab can establish or verify the accuracy of the measurement response.

AMR Verification

- AMR verification can be based on demonstrating linearity either by measuring recoveries or by demonstrating expected proportions


Materials for AMR Verification

- Linearity material of appropriate matrix
- PT (or PT validated) materials
- Previously tested clinical specimens
- Previously tested clinical specimens, altered by admixture, dilution, spiking
- Reference materials with matrix and target values appropriate for the method
- Method product calibrator(s)

Analytical Measurement Range Verification

CONCENTRATION EXPECTED

Based on: (a) admixture ratio (b) assigned value

AMR Verification



1. Specimens (can use patients) with low and high values
2. Mix 1:1 to get midpoint value
3. Mix low:mid and mid:high to get additional values
4. Expected values are:

<u>Sample</u>	<u>Proportion</u>	<u>Concentration expected</u>
Low	0%	20 mg/dL
0.75 L + 0.25 H	25%	160 mg/dL
0.50 L + 0.50 H	50%	299 mg/dL
0.25 L + 0.75 H	75%	438 mg/dL
High	100%	578 mg/dL

Criteria for acceptability = recovery of concentration in proportion to the amount present (linear relationship)

AMR Verification

Linear relationship verifies recovery of correct concentration relationship over AMR.

Clinically Reportable Range

- Range of analyte values that will be reported as a quantitative result, allowing for specimen dilution or other pretreatment to extend the AMR
- Policy based on the laboratory director's medical judgment regarding the clinical use requirements for a laboratory test

Clinically Reportable Range

Example: hCG may have a CRR with both low and high limits.

- AMR is 3-1,000 mIU/mL (method specification).
- CRR is 5-1,000,000 mIU/mL (clinical usefulness).
- Results of <3, 3, or 4 are reported as "<5 mIU/mL".
- Results >1,000 are diluted and rerun to obtain quantitative values up to 1,000,000 mIU/mL.
- Results >1,000,000 are reported as ">1,000,000 mIU/mL".
- NOTE: CRR limits are the judgment of the lab director and may be different for different labs

Reference Limits

- Not required for laboratory to establish its own reference limits, but must verify that limits it uses are appropriate for patient population served
- Should review how any published limits were established before considering whether to try to adopt them

Adopted Reference Limits

- Manufacturer suggested
- Reference laboratory
- Published articles
- Neighboring laboratory
- Previous reference limits in own lab

Reference Limits

- Validation protocol outlined in CLSI C28-A2
- Select 20 representative healthy individuals and do test; if ≤ 2 outside proposed limits, validated
- If > 2 outside, can repeat with another 20, and accept if ≤ 2 outside (not worth repeat if > 5 outside proposed limits)

Example

- Old CI reference limits 98-104; data below from 20 patients. What would you do?

105	103	107	102	106
101	103	107	101	100
102	108	104	105	103
100	107	102	105	101

Reference Limits

- If not validated, will need to establish own reference limits
- Briefly, need 120 individuals (200 if does not follow gaussian distribution) to be confident in accuracy
- Need to have criteria for exclusion (may differ for certain tests)


Resources

- CAP Laboratory General Checklist
- CLSI EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods (2004)
- CLSI EP6-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach (2003)

Resources

- CAP Laboratory General Checklist
- CLSI EP5-A2: Evaluation of Precision Performance of Quantitative Measurement Methods (2004)
- CLSI EP6-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach (2003)

Sensitivity measurement

Method Sensitivity

- The *analytical* sensitivity of a method refers to the lowest concentration of analyte that can be reliably detected.
- The most common definition of sensitivity is the analyte concentration that will result in a signal two or three standard deviations above background.

Clinical Sensitivity Vs Analytical Sensitivity

Other measures of sensitivity

- Limit of Detection (LOD) is sometimes defined as the concentration producing an $S/N > 3$.
 - In drug testing, LOD is customarily defined as the lowest concentration that meets all identification criteria.
- Limit of Quantitation (LOQ) is sometimes defined as the concentration producing an $S/N > 5$.
 - In drug testing, LOQ is customarily defined as the lowest concentration that can be measured within $\pm 20\%$.

Limit of Detection

- The limit of detection is the point at which a measured value is larger than the uncertainty associated with it. **It is the lowest concentration of analyte in a sample that can be detected but not necessarily quantified.** The limit of detection is frequently confused with the sensitivity of the **method**. The sensitivity of an analytical method is the capability of the method to discriminate small differences in concentration or mass of the test analyte. In practical terms, sensitivity is the slope of the calibration curve that is obtained by plotting the response against the analyte concentration or mass.

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- The limit of quantitation is the **minimum injected amount that produces quantitative measurements in the target matrix with acceptable precision** in chromatography, typically requiring peak heights 10 to 20 times higher than the baseline noise.
- . A number of samples with decreasing amounts of the analyte are injected six times. The calculated RSD percent of the precision is plotted against the analyte amount. The amount that corresponds to the previously defined required precision is equal to the limit of quantitation. It is important to use not only pure standards for this test but also spiked matrices that closely represent the unknown samples.

Verify Analytic Sensitivity

- Also called lower detection limit
- Involves two steps: determination of values obtained with blank samples, and values obtained with low level positive samples
- Blank samples often use the zero calibrator for an assay, as it can be difficult to find samples truly missing the substance of interest

Verify Analytic Sensitivity

- Low level positive samples identified at or only slightly above the manufacturer's stated analytical sensitivity; for uncommon tests, may need to contact another laboratory doing the test for help in finding low positive samples
- CLSI EP17-A describes actual procedure
 - Run 20 blanks; if < 3 exceed stated blank value, accept that value
 - Run 20 low patient samples near the detection limit; if at least 17 are above the blank value, the detection limit is verified

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Evaluating Clinical Performance of laboratory tests

- The *clinical performance* of a laboratory test defines how well it predicts disease
- The *sensitivity* of a test indicates the likelihood that it will be positive when disease is present
- The *specificity* of a test indicates the likelihood that it will be negative when disease is absent

Clinical Sensitivity

If TP as the number of “true positives”, and FP is the number of “false positives”, the *sensitivity* is defined as:

$$Sensitivity = \frac{TP}{TP + FN} \cdot 100$$

Example

Of 25 admitted cocaine abusers, 23 tested positive for urinary benzoylecgonine and 2 tested negative.

What is the sensitivity of the urine screen?

$$\frac{23}{23+2} \cdot 100 = 92\%$$

Clinical Specificity

If TN is the number of “true negative” results, and FP is the number of falsely positive results, the *specificity* is defined as:

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